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<https://www.mathworks.com/help/daq/acquire-data-in-the-background.html>

https://it.mathworks.com/help/matlab/matlab_oop/learning-to-use-events-and-listeners.html?searchHighlight=Events%20and%20Listeners%20%E2%80%94%20Concepts&s_tid=doc_srchtile

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Esempi_background_DAQ

a. Acquire Dati in Background

A background acquisition depends on events and listeners to allow your code to access data as the hardware acquires it and to react to any errors as they occur. For more information, see Events and Listeners — Concepts in the MATLAB Object-Oriented Programming documentation.

Use events to acquire data in the background. In this example, you acquire data from an NI 9205 device with ID cDAQ1Mod1 using a listener and a DataAvailable event.

Listeners execute a callback function when notified that the event has occurred. Use `Session.addListener` to create a listener object that executes your callback function.

Create an NI session object and an analog input 'Voltage' channel on cDAQ1Mod1:

```
s = daq.createSession('ni');  
addAnalogInputChannel(s, 'cDAQ1Mod1', 'ai0', 'Voltage');
```

Add the listener for the DataAvailable event and assign it to the variable lh:

```
lh = addlistener(s, 'DataAvailable', @plotData);
```

For more information on events, see Events and Listeners — Concepts in the MATLAB Object-Oriented Programming documentation.

Create a simple callback function to plot the acquired data and save it as `plotData.m` in your working directory:

```
function plotData(src,event)  
    plot(event.TimeStamps, event.Data)  
end
```

Here, `src` is the session object for the listener and `event` is a `daq.DataAvailableInfo` object containing the data and associated timing information.

Acquire the data and see the plot update while MATLAB® is running:

```
startBackground(s);
```

When the operation is complete, delete the listener:

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```
Type('stopWhenExceedOneV.m')
```

```
function stopWhenExceedOneV(src, event)
    if any(event.Data > 1.0)
        disp('Event listener: Detected voltage exceeds 1, stopping acquisition')
        % Continuous acquisitions need to be stopped explicitly.
        src.stop()
        plot(event.TimeStamps, event.Data)
    else
        disp('Event listener: Continuing to acquire')
    end
end
```

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```
while s.IsRunning
    pause(0.5)
    fprintf('While loop: Scans acquired = %d\n', s.ScansAcquired)
end
```

```
fprintf('Acquisition has terminated with %d scans acquired\n', s.ScansAcquired);
```

```
Event listener: Continuing to acquire
While loop: Scans acquired = 1007
Event listener: Detected voltage exceeds 1, stopping acquisition
While loop: Scans acquired = 1207
Acquisition has terminated with 1207 scans acquired
```